

Semi-Flex™ Elbow/Wrist Coil

ScanMed's Semi-Flex™ Elbow/Wrist Coil is fully optimized to provide high quality soft tissue and bone imaging of the extremities. The unique elliptical shape combined with our proprietary Semi-Flex™ technology brings the array elements closer to the targetted regions for precise imaging, higher SNR, and more uniformity.

The semi-flexible foam housing and 72 inch cable allows the technologist to easily position the coil on the extremities, enhancing comfort and allowing for better patient compliance and faster throughput.

This versatile coil images the following anatomies:

- Elbow
- Wrist
- Forearm
- Hand
- Fingers
- Forefoot
- Toes



SPECIFICATIONS:

- Field Strength: 1.5T, 3T
- 8 or 16 elements for 8 or 16 channel systems
- Cable length: 183 cm (72 in)
- Weight: 1.4 kg (3 lb)

Includes foam positioning aid



*Images of a wrist and elbow using our
8-Channel Semi-Flex™ Elbow/Wrist Coil*

**Contact us for
more information**



✉ sales@scanmed.com

🌐 www.scanmed.com

☎ 402.934.2650

📄 402.778.9699