

# Knee/Foot Coil

ScanMed's Knee/Foot Coil has a unique chimney design that allows for versatile 360° imaging of the knee, foot, ankle, and distal extremities without repositioning. The highly-sensitive, optimized elements of the provide broad coverage, improved SNR, and impressive uniformity for all targeted anatomies.

The Knee/Foot Coil is the only coil on the market to come with a wide bore separable housing to accommodate larger patients without compromising SNR or uniformity. The coil set also includes a soft foam rest for both legs facilitating patient comfort and easy image acquisition.



## SPECIFICATIONS:

- Field Strength: 1.5T, 3T
- 8 elements for 8 or 16 channel systems
- Cable length: 127 cm (48 in)
- Weight: 5.5 kg (12)



*Images of the knee and foot using our 8-Channel Knee/Foot Coil*

Contact us for  
more information



✉ sales@scanmed.com

🌐 www.scanmed.com

☎ 402.934.2650

📄 402.778.9699