

Pediatric Blanket™ Coil

ScanMed's Pediatric Blanket™ Coil is the first-of-a-kind MRI coil in a blanket. This innovative coil was designed to accommodate situations where skin sensitivity, minimal pressure, and patient comfort are of vital importance. The Pediatric Blanket™ Coil perfectly follows the contours of a healthy or tortuous spine and can safely image the full range of pediatric anatomies.


Inside the lightweight, quilt-like enclosure is a highly flexible MRI antenna set that wraps gently around a pediatric body, offering similar comforts to swaddling. A soft antimicrobial fleece cover provides another layer of coziness, keeping the infant fully relaxed and still during scan time.

The Pediatric Blanket™ Coil is also extremely versatile and can be easily used for scanning a variety of adult anatomies.



SPECIFICATIONS:

- Field Strength: 1.5T, 3T
- 8, 16, or 32 elements for 8, 16, or 32 channel systems
- Cable length: 150 cm (59 in)
- Weight: 1.3 kg (3



Pediatric images provided by our 16 Channel Pediatric Blanket™ Coil

Contact us for more information



✉ sales@scanmed.com

🌐 www.scanmed.com

☎ 402.934.2650

📄 402.778.9699