

Wrist Coil

ScanMed's dedicated Wrist Coil comes in two element configurations. In the 8-element Wrist Coil, all 8 antennae surround the wrist and are optimized in length to cover only the bone and joint region. In the 16-element Wrist Coil, there are 8 additional antennae added in an overlapping design to fully surround the hands through the first joints of the fingers providing homogenous imaging volume. Both configurations have demonstrated significant image quality improvement in a FOV as small as 8cm.

The Wrist Coil comes with two separable bases for positioning either by the side of the patient or above their head when in the superman pose. With versatile positioning options and a soft foam interior lining, the Wrist Coil also offers incredible patient comfort for faster and more precise image acquisitions.



SPECIFICATIONS:

- Field Strength: 1.5T, 3T
- 8 or 16 elements for 8 or 16 channel systems
- Cable length: 109 cm (43 in)
- Weight: 3.1 kg (7 lb)



Images of the wrist using our 8-Channel dedicated Wrist Coil

Contact us for
more information



✉ sales@scanmed.com

🌐 www.scanmed.com

☎ 402.934.2650

📄 402.778.9699